

Wednesday 19 June 2019 - Morning

A Level Chemistry B (Salters)

H433/03 Practical skills in chemistry

Practical Insert

Time allowed: 1 hour 30 minutes



INSTRUCTIONS

• Do not send this insert for marking; it should be retained in the centre or destroyed.

INFORMATION

• This document consists of 4 pages. Any blank pages are indicated.

Iron in spinach

Spinach has often been regarded as an excellent source of dietary iron.

Below a student describes an investigation to determine the mass of iron contained in a typical portion of spinach used in a meal.

Introduction

The amount of iron, as Fe^{2+} , in spinach can be found by titration with potassium manganate(VII) solution.

Manganate(VII), MnO_4^- , is a strong oxidising agent. It accepts electrons easily, and is reduced to colourless manganese(II) ions according to the half-equation below:

$$\mathrm{MnO_4}^-(\mathrm{aq})$$
 + 8H⁺(aq) + 5e⁻ \rightarrow $\mathrm{Mn^{2+}(aq)}$ + 4H₂O(I) (purple) (colourless)

The electrons are provided by reducing agents such as iron(II) salts:

$$Fe^{2+}(aq) \rightarrow Fe^{3+}(aq) + e^{-}$$

As a result, manganate(VII) can be used in acidic solution to determine the number of moles of reducing agent, e.g. Fe²⁺, present.

Manganate(VII) is added from a burette to a solution of Fe²⁺ ions and is decolourised immediately. As soon as the Fe²⁺ ions are used up, the next drop of manganate(VII) is not decolourised, and so the solution in the conical flask goes pale pink. The end-point of the titration is the first permanent appearance of this pale pink colour. Manganate(VII) is therefore self-indicating and no other indicator is needed.

The acid used to provide $H^+(aq)$ is dilute sulfuric acid; this should always be in excess or else insoluble brown MnO_2 will form.

Getting the Fe²⁺ ions into solution

Approximately 5g portions of spinach were immersed in dilute sulfuric acid for various amounts of time. The solutions were filtered and 25 cm³ portions were titrated with the standard potassium manganate(VII) solution.

© OCR 2019 H433/03 Jun19

Method

- 1. Four samples of approximately 5g of the spinach leaves provided were weighed by difference, accurately, using a 2 decimal place balance. All the weighings were recorded.
- 2. Each weighed sample of spinach was added to about 100 cm³ of sulfuric acid in a beaker and allowed to stand for various amounts of time.
 - After standing each sample was filtered into a 250 cm³ volumetric flask. The original beakers were washed several times with de-ionised water and the washings transferred to the flask. The solution was made up to the mark with de-ionised water.
- 3. 25 cm³ of one of the solutions was pipetted into a conical flask.
- 4. The above solution was titrated against a $5.0 \times 10^{-6} \, \text{mol dm}^{-3}$ solution of KMnO₄ from a burette until at least two concordant results were obtained.
- 5. Steps 3, 4 and 5 were repeated with each of the sample solutions.

Results and Analysis

Weighings

	Mass of weighing boat/g	Mass of spinach + weighing boat/g	Mass of spinach/g
Sample 1	1.43	6.75	5.32
Sample 2	1.43	6.98	5.55
Sample 3	1.43	6.40	4.97
Sample 4	1.43	6.53	5.10

Table 1

Titrations

		Sample 1	Sample 2	Sample 3	Sample 4
	Time/mins	30	60	90	120
Rough titre	Initial vol/cm ³	0.00	0.00	0.00	0.00
	Final vol/cm ³	6.80	7.25	7.70	8.20
	Titre/cm ³	6.80	7.25	7.70	8.20
Titre 1	Initial vol/cm ³	7.00	8.00	8.00	10.00
	Final vol/cm ³	13.80	15.10	15.55	18.05
	Titre/cm ³	6.80	7.10	7.55	8.05
Titre 2	Initial vol/cm ³	15.00	16.00	16.00	20.00
	Final vol/cm ³	21.75	23.15	23.50	27.95
	Titre/cm ³	6.75	7.15	7.50	7.95



Copyright Information

OCR is committed to seeking permission to reproduce all third-party content that it uses in its assessment materials. OCR has attempted to identify and contact all copyright holders whose work is used in this paper. To avoid the issue of disclosure of answer-related information to candidates, all copyright acknowledgements are reproduced in the OCR Copyright Acknowledgements Booklet. This is produced for each series of examinations and is freely available to download from our public website (www.ocr.org.uk) after the live examination series.

If OCR has unwittingly failed to correctly acknowledge or clear any third-party content in this assessment material, OCR will be happy to correct its mistake at the earliest possible opportunity.

For queries or further information please contact The OCR Copyright Team, The Triangle Building, Shaftesbury Road, Cambridge CB2 8EA.

OCR is part of the Cambridge Assessment Group; Cambridge Assessment is the brand name of University of Cambridge Local Examinations Syndicate (UCLES), which is itself a department of the University of Cambridge.

© OCR 2019 H433/03 Jun19